

A new idea that really holds water (and soap)

The Vernaclean Bowl is unique, the only maceratable pulp bowl that holds warm water and detergent.

Ideal for: Bed baths | Hand washing | General surface cleaning | Wound care | Shaving | Contenance care

Now available
from NHS
supplies

Prevents infection outbreaks
associated with plastic bowls

Saves nurse time and money



Maximum capacity 4 litres

Fill level 2 litres

Holds warm water and detergent for up to 2 hours

Suitable for use with detergents:

Including: Soap | Shampoo | Conditioner | Washing up liquid

User feedback from evaluations involving 92 nurses*:

98.5% agreed that the Vernaclean bowl will reduce the likelihood of infections.

93.2% agreed that the product saved nurse time.

Patients liked it as they got a clean bowl every time.

* Evaluations across four NHS trusts including: Trafford General Hospital, Glasgow Royal Infirmary, Leeds Teaching Hospital, Cirencester Hospital.

Now available from NHS supplies or direct from Vernacare

Order Codes:

Description:	NHS Catalogue Code (England)	Vernacare Code (Scotland & Wales)
Case of 50 Vernaclean bowls	FSL307	118AA050
Vernaclean Rack	-	510ZS001

Instructions for use:



Place upright on top of blades to macerate

1 Remove bowl from rack

2 Fill with warm water (max temp 45°C). Add detergent

3 Dispose via your macerator:
Vernacare: up to 2 at a time
Haigh: up to 1 at a time